


Introduction

Iron oxide nanoparticles (NPs) have been used for a variety of *in vivo* and *ex vivo* applications in the biomedical sciences. Moreover, when intended for *in vivo* clinical applications, NPs need to meet rigorous requirements to ensure safety as well as bio-functionality, including blood circulation time and specificity for cellular targets. PrecisionMRX[®] NPs are extensively characterized superparamagnetic NPs composed of 25 nm magnetite cores that are currently used in a variety of *in vivo* applications including non-invasive *in vivo* diagnosis of cancer, Magnetic Particle Imaging, MRI, and magnetic hyperthermia.

Objective

Conduct development and pre-clinical functionality studies of anti-HER2 antibody (mAb) conjugated NPs for *in vitro* and *ex vivo* detection of HER2+ tumor cells by Magnetic Relaxometry (MRX).

Superparamagnetic Relaxometry (SPMR) Technology




SPMR only detects NPs bound to cells/tissues

1. Inject NPs
2. Apply small magnetizing pulse
3. Turn off field
4. NPs relax to their equilibrium states.

- **Brownian** motion of unbound NPs (fast)
- **Néel** relaxation of NPs bound to cells (slow and measurable)

PrecisionMRX[™] NANOPARTICLES




Methods

Anti-HER2 NP functionalization and characterization

PrecisionMRX[®] NPs were encapsulated by a layer of polymer and then functionalized with carboxylate (COO⁻) surface. PEG + anti-HER2 mAb were subsequently conjugated onto the polymer surface. Size of resulting NPs were measured by DLS. Bound and free mAb were determined via ELISA.

Surface	Diameter	PDI	# of Ab/NP	% of free Ab
PEG + anti-HER2	70-80 nm	<0.10	3-5	<10%




Specific binding *in vitro*

- A variety of cell lines with different levels of HER2 expression were incubated with 100ug of anti-HER2 mAb NPs overnight.
- Cells were washed, harvested, centrifuged, and pellets were subsequently measured on the MRX instrument.
- Cell competition study was done by pre-incubate cells with free anti-HER2 antibody.


Specific binding *in vivo*

- Anti-HER2 mAb conjugated NPs or PEGylated NPs were injected into BT474 (HER2 (3+)) and MCF-7 (HER2 (1/0+)) tumor bearing mice by intraperitoneal (IP) or peritumoral (PT) delivery.
- For *in vivo* competition, free anti-HER2 mAb were injected 24hr prior to NP injection.
- After 24 hr post NP injection, mice were euthanized and tumors and organs were excised for *ex vivo* MRX measurement




Results – In Vitro Specificity, Sensitivity

Cell Competition Assay




Cell Titration Assay




Nanoparticles can clearly distinguish high, medium and low HER2 expressing cell lines, demonstrating specificity and sensitivity

Results – In Vivo Specificity

Tumored mice: MRX signals of resected tumors




Non-Tumored mice: MRX signals of major organs




- MRX tumor binding signal can be competed out by free anti-HER2 mAb indicating specificity
- Anti-HER2 mAb conjugated NP generated higher binding signal than PEGylated NP indicating specificity
- All major organs except liver (clearance) have no significant NP accumulation indicating anti-HER2 mAb NP is safe

Dual Flank Tumors of BT474 (HER2 3+) and MCF7 (HER2 1/0+)

MRX signals of resected tumors



IHC of Tumor (Anti-HER2 Secondary Ab)



Dual flank tumor study demonstrated that BT474 tumor generated much higher binding signals compared to MCF7 tumor. These results were confirmed by the presence/absence of anti-HER2 mAb in BT474 and MCF7 tumor respectively in IHC study.

Conclusions and Future Work

Together, these results suggest that our anti-HER2 antibody conjugated nanoparticles are safe; can provide targeted and specific delivery to cancerous tissue *in vivo* and generate measurable signal on our MRX detection instrument. These studies lay out ground work for our future human clinical study for breast cancer detection.